

Morinaga Food Allergen Lateral Flow  
Rapid Test Pro II for Crustacean

Rev. Mar. 2024

Lot number from #PR001

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result	
<b>Cereals</b>	Rice powder	-	<b>Seeds</b>	Peanut	-	<b>Fruits</b>	Orange	-	
	Polished rice	-		Roasted almond	-		Kiwi	-	
	Sprouted brown rice	-		Roasted cashew nut	-		Apple	-	
	Glutinous rice	-		Roasted macadamia nut	-		Peach	-	
	Buckwheat	-		Roasted pistachio nut	-		Banana	-	
	Japanese barnyard millet	-		Roasted walnut	-	Coconut milk	-		
	Foxtail millet	-		Roasted pecan nut	-	Coffee bean	-		
	Common millet	-		White sesame	-	Cocoa powder	-		
	Wheat	-		Black sesame	-	Green tea	-		
	Rye	-		Salmon	-	Roasted laver	-		
	Barley	-	Mackerel	-	Hijiki (boiled)	-			
	Malt	-	Squid	-	Brown seaweed (raw)	-			
	Oat	-	Octopus	-	Kelp 200-fold extract	-			
	Adlay	- *1	<b>Sea foods</b>	Black tiger prawn	+ *2	<b>Spices</b>	Cumin	-	
	Corn	-	Snow crab	+	Poppy seeds		-		
Corn flour	-	Scallop adductor	-	Black pepper	-				
Amaranth	-	Pollack roe	-	White pepper	-				
Quinoa	-	Salmon roe	-	Red pepper	-				
<b>Tubers· Starch</b>	Sago palm powder	-	<b>Meats</b>	Beef	-		Ginger	-	
	Yam	-		Chicken	-		Cinnamon	-	
	Potato starch	-		Pork	-		Coriander	-	
<b>Beans</b>	Soybean	-	<b>Eggs</b>	Chicken egg	-		<b>Polysaccharide Thickener</b>	Turmeric	-
	Red bean (azuki)	-		Processed chicken egg	-			Low methoxyl pectin 50-fold extract	-
	Common bean	-	<b>Dairy products</b>	Dried chicken egg	-	<b>Others</b>	Guar gum 1000-fold extract	-	
	Concord pole	-		Milk	-		Gelatin	-	
	Broad bean	-		Skim milk	-				
		<b>Green stuff</b>	Tomato	-					
			Onion	-					
		<b>Fungus</b>	Spinach	-					
			Mushroom	-					

Test Result	
-	Negative
+	Positive

\*1 The Test Solution may be difficult to flow on a test stick.

\*2 The test result may be negative if it is included in high concentrations.

These samples were tested by following the instruction of Rapid Test Pro II (using Recommended Extraction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.